

# IINTS-AF SDK Technical Reference Manual

Research Use Only — Not for Clinical Care

Generated: 2026-03-01 12:59 UTC

## IINTS-AF SDK Technical Reference Manual

Version 0.1.19 | Python SDK

**\*\*PRE-CLINICAL USE ONLY - NOT FOR PATIENT CARE\*\***

This SDK is intended for research, simulation, and algorithm validation.  
It has NOT received FDA clearance or CE marking for clinical use.

### How to Use This Manual (Read First)

This manual is long. Use this map to find what you need fast:

**\*\*Full SDK overview:\*\*** ``docs/COMPREHENSIVE_GUIDE.md``

**\*\*CLI reference:\*\*** ``docs/TECHNICAL_README.md``

**\*\*Research track (predictor training):\*\*** ``research/README.md``

**\*\*Notebooks:\*\*** ``examples/notebooks/README.md``

Recommended notebook order (best for new users):

``00_Quickstart.ipynb`` — run a full simulation

``01_Presets_and_Scenarios.ipynb`` — scenarios + presets

``02_Safety_and_Supervisor.ipynb`` — safety checks

``03_Audit_Trail_and_Report.ipynb`` — audit trail + PDF

``04_Baseline_and_Metrics.ipynb`` — metrics & baselines

``05_Devices_and_HumanInLoop.ipynb`` — pumps/sensors + manual override

``06_Optional_Torch_LSTM.ipynb`` — predictor training

``07_Ablation_Supervisor.ipynb`` — safety ablation

``08_Data_Registry_and_Import.ipynb`` — data import + registry

Quick task routing:

**\*\*Run a simulation fast:\*\*** Section 2.2

**\*\*Customize safety limits:\*\*** Section 4.2

**\*\*Generate audit report:\*\*** Section 5.4

**\*\*Train an AI predictor:\*\*** Section 8.2 + ``research/README.md``

### Table of Contents

#### 1. Executive Summary

Overview and key capabilities

Intended use and audiences

Safety-first philosophy

#### 2. Getting Started (Step-by-Step)

##### 2.1 Installation Guide

##### 2.2 Your First Simulation (60 seconds)

## 2.3 Understanding the Output Files

## 2.4 Creating Custom Algorithms

## 2.5 Adding Stress Tests

## 2.6 Benchmarking Against Baselines

## 2.7 Using Preset Scenarios

## 2.8 Importing Real CGM Data

## 2.9 Generating Custom Reports

## 2.10 Next Steps

## 3. Architecture Overview

### 3.1 System Components

### 3.2 Data Flow Diagram

### 3.3 Safety Layer Integration

## 4. Safety Architecture (Critical)

### 4.1 Design Philosophy

### 4.2 SafetyConfig Configuration

### 4.3 9 Safety Checks Explained

### 4.4 Safety Levels and Interventions

### 4.5 Input Validation Rules

### 4.6 Simulation Termination Conditions

## 5. Tutorials and Cookbook

### 5.1 24-Hour Simulation Walkthrough

### 5.2 Building an ML-Hybrid Algorithm

### 5.3 Batch Experiment Execution

### 5.4 Audit Trail Analysis

### 5.5 Custom Safety Thresholds

### 5.6 Pump Emulator Benchmarking

### 5.7 Live Streaming Simulation

### 5.8 Reproducible Runs for Publications

## 6. API Reference

### 6.1 Core Classes and Interfaces

### 6.2 Algorithm Development Guide

### 6.3 Simulator Configuration

### 6.4 Patient Profile Customization

### 6.5 Device Models (Sensor/Pump)

### 6.6 High-Level API Functions

## 7. Practical Examples

### 7.1 Complete Algorithm Example

### 7.2 CLI vs Python API Comparison

### 7.3 Stress Testing Patterns

### 7.4 Human-in-the-Loop Integration

### 7.5 Data Import Workflows

## 8. Advanced Topics

### 8.1 Commercial Pump Emulators

### 8.2 Dataset Registry Usage

### 8.3 Reproducibility Techniques

### 8.4 Performance Profiling

### 8.5 Custom Metrics Calculation

## 9. Troubleshooting

### 9.1 Common Installation Issues

### 9.2 Simulation Problems

### 9.3 Algorithm Development Tips

### 9.4 Data Import Solutions

### 9.5 Performance Optimization

## 10. Quick Reference

### 10.1 Essential CLI Commands

### 10.2 Python Code Snippets

### 10.3 Safety Thresholds Cheatsheet

### 10.4 Clinical Metrics Targets

## 11. Glossary of Terms

## 1. Executive Summary

The IINTS-AF SDK (Intelligent Insulin Titration System for Artificial Pancreas) is a **\*\*safety-first simulation and validation platform\*\*** for insulin dosing algorithms targeting closed-loop insulin delivery research.

### Key Capabilities

[OK] **\*\*Plug-and-play algorithm architecture\*\*** - Implement one method, get full simulation

[OK] **\*\*9-layer Independent Safety Supervisor\*\*** - Deterministic override guarantees

[OK] **\*\*Realistic device models\*\*** - CGM sensor and insulin pump error simulation

[OK] **\*\*Commercial pump emulators\*\*** - Medtronic 780G, Omnipod 5, Tandem Control-IQ

[OK] **\*\*Clinical metrics\*\*** - TIR, GMI, CV, LBGI, HBGI per ATTD/ADA guidelines

[OK] **\*\*Complete audit trail\*\*** - JSONL + CSV + JSON summary with integrity hashing

[OK] **\*\*PDF clinical reports\*\*** - Visual reports with glucose traces and safety summaries

[OK] **\*\*Benchmark mode\*\*** - Head-to-head algorithm comparison

[OK] **\*\*Real-world data import\*\*** - Dexcom, Libre, Nightscout, AIDE/PEDAP, AZT1D, HUPA-UCM

[OK] **\*\*Optional AI predictor\*\*** - Proactive glucose forecasting

[OK] **\*\*Reproducible runs\*\*** - Seeded randomness and signable manifests

### Intended Use

This SDK is intended for:

**\*\*Pre-clinical algorithm validation\*\***

**\*\*Academic research\*\***

**\*\*Educational purposes\*\***

**\*\*Regulatory submission preparation\*\***

**\*\*NOT intended for:\*\***

- Direct patient care
- Clinical decision-making
- Medical device deployment without regulatory review

## 2. Getting Started

### 2.1 Installation Guide

#### Option 1: Install from PyPI (Recommended)

```
# Create project folder
mkdir my-aps-research && cd my-aps-research

# Create virtual environment (recommended)
python3 -m venv .venv
source .venv/bin/activate # macOS/Linux
# .venv\Scripts\activate # Windows

# Install SDK
pip install iints-sdk-python35

# Verify installation
iints --help
```

#### Option 2: Development Install

```
git clone https://github.com/python35/IINTS-SDK.git
cd IINTS-SDK
pip install -e ".[dev]"
```

#### Option 3: With Research Extras (AI Predictor)

```
pip install iints-sdk-python35[research]
```

### 2.2 Your First Simulation (60 seconds)

#### Using CLI (Quickstart)

```
# Create quickstart project
iints quickstart --project-name iints_quickstart
cd iints_quickstart

# Run clinic-safe preset
iints presets run --name baseline_t1d --algo algorithms/example_algorithm.py
```

#### Using Python API

```
from iints import run_simulation
from iints.core.algorithms.pid_controller import PIDController

# Run 12-hour simulation
results = run_simulation(
    algorithm=PIDController(),
    duration_minutes=720,
    seed=42
)
```

```
print(f"Results saved to: {results['results_csv']}")
print(f"Report generated: {results['clinical_report']}")
```

**\*\*What this does:\*\***

- Creates a virtual patient with default parameters
- Runs PID controller algorithm for 12 hours
- Generates results CSV, clinical report PDF, and audit trail
- Compares against baseline algorithms automatically

**\*\*Sanity check (first run):\*\***

- `results/results.csv` exists and grows during the run
- `results/clinical\_reports/` contains a PDF
- Console shows no safety contract violations

**\*\*Expected runtime:\*\***

- Laptop CPU: ~30-90 seconds for the quickstart preset
- GPU not required

## 2.3 Understanding the Output Files

After running a simulation, you'll find these files in the `results/` folder:

File	Description
----- -----	
`results.csv`	Every simulation step with glucose, insulin, IOB, COB, safety events
`clinical_report.pdf`	Visual report with charts and clinical metrics
`config.json`	Exact configuration used (for reproducibility)
`run_metadata.json`	Run ID, seed, platform info, timestamps
`run_manifest.json`	SHA-256 hashes of all files (integrity verification)
`audit/audit_trail.csv`	Detailed per-step audit trail
`audit/safety_summary.json`	Safety interventions summary
`baseline/pid_results.csv`	PID controller baseline comparison
`baseline/standard_pump.csv`	Standard pump baseline comparison

**\*\*Quick interpretation:\*\***

- Start with `clinical\_report.pdf` for a human-readable summary.
- Use `results.csv` for plotting and deeper analysis.
- Use `audit/safety\_summary.json` to explain *\*why\** the supervisor intervened.

**\*\*Data consistency checks:\*\***

- `glucose\_actual\_mgdl` should be in 40-400 mg/dL range.
- `patient\_iob\_units` and `patient\_cob\_grams` should be  $\geq 0$ .
- Large jumps ( $>60$  mg/dL in 5 min) usually indicate data issues.

**\*\*Example: Loading Results in Python\*\***

```
import pandas as pd

# Load main results
df = pd.read_csv("results/your_run/results.csv")
print(df.head())
```

```
# Load safety summary
import json
with open("results/your_run/audit/safety_summary.json") as f:
    safety = json.load(f)
print(f"Total interventions: {safety['intervention_count']}")
```

## 2.4 Creating Your First Custom Algorithm

### Step 1: Generate Template

```
iints new-algo --name MyAlgorithm --output-dir algorithms/
```

### Step 2: Implement Logic

```
# algorithms/my_algorithm.py
from iints.api.base_algorithm import InsulinAlgorithm, AlgorithmInput

class MyAlgorithm(InsulinAlgorithm):
    def get_algorithm_metadata(self):
        return {
            "name": "My Custom Algorithm",
            "version": "0.1.0",
            "description": "My first insulin dosing algorithm"
        }

    def predict_insulin(self, data: AlgorithmInput) -> dict:
        # Simple example: basal rate + correction
        basal = 0.9 # U/hr
        correction = 0.0

        # Add correction if glucose is high
        if data.current_glucose > 180:
            correction = (data.current_glucose - 180) / 50 # ISF of 50

        return {
            "basal_insulin": basal,
            "bolus_insulin": correction,
            "reason": f"Basal {basal}U + correction {correction}U for glucose {data.current_glucose}mg/dL"
        }
```

### Step 3: Test Your Algorithm

```
# Run with your custom algorithm
iints run --algo algorithms/my_algorithm.py --duration 1440
```

## 2.5 Adding Stress Tests

Add realistic challenges to test algorithm robustness:

```
from iints.core.simulator import Simulator, StressEvent

sim = Simulator(algorithm=MyAlgorithm(), patient_config="default")

# Add meal at 8:00 AM (480 minutes)
sim.add_stress_event(StressEvent(
    start_time=480,
```

```

    event_type="meal",
    value=60 # 60g carbs
))

# Add exercise at 6:00 PM (1080 minutes)
sim.add_stress_event(StressEvent(
    start_time=1080,
    event_type="exercise",
    value=30 # 30 minutes moderate exercise
))

# Add sensor noise
sim.add_stress_event(StressEvent(
    start_time=300,
    event_type="sensor_noise",
    value=15.0 # 15 mg/dL standard deviation
))

results = sim.run_batch(duration_minutes=1440)

```

#### **\*\*Available Stress Events:\*\***

- `meal`: Carbohydrate intake (value = grams)
- `exercise`: Physical activity (value = minutes)
- `sensor\_noise`: CGM noise (value = std deviation)
- `sensor\_dropout`: CGM signal loss (value = probability)
- `pump\_failure`: Insulin delivery failure (value = probability)
- `hormonal\_change`: Dawn phenomenon simulation

## 2.6 Benchmarking Against Baselines

Compare your algorithm against standard approaches:

```

# CLI benchmark
iints benchmark \
  --algo-to-benchmark algorithms/my_algorithm.py \
  --patient-configs-dir src/iints/data/virtual_patients \
  --scenarios-dir scenarios \
  --output-dir results/benchmark

# Python API benchmark
from iints.analysis.benchmark import run_benchmark

results = run_benchmark(
    algorithm_path="algorithms/my_algorithm.py",
    patient_configs=["default", "adolescent", "insulin_resistant"],
    scenarios=["baseline", "meal_challenge", "exercise_stress"],
    duration_minutes=1440
)

```

## 2.7 Using Preset Scenarios

Clinic-safe scenarios for reproducible testing:

```

# List available presets
iints presets list

```

```
# Run a specific preset
iints presets run --name hypo_prone_night --algo algorithms/my_algorithm.py
```

**\*\*Available Presets:\*\***

- `baseline\_t1d`: Standard Type 1 diabetes profile
- `hypo\_prone\_night`: Nighttime hypoglycemia risk
- `hyper\_challenge`: Post-meal hyperglycemia
- `pizza\_paradox`: Delayed carb absorption
- `midnight\_crash`: Nocturnal hypoglycemia
- `exercise\_stress`: Physical activity impact
- `sensor\_noise`: CGM accuracy challenges

## 2.8 Importing Real CGM Data

Test against real patient data:

```
# From CSV file
iints import-data --input-csv my_cgm_export.csv \
  --data-format dexcom \
  --output-dir results/imported

# Python API
from iints.data.importer import import_cgm_csv
from iints.core.simulator import Simulator

result = import_cgm_csv(
    "my_cgm_export.csv",
    data_format="dexcom", # or "libre", "generic"
    scenario_name="Patient A - Week 1",
)

sim = Simulator(
    algorithm=MyAlgorithm(),
    scenario=result.scenario,
)
```

**\*\*Supported Formats:\*\***

- Dexcom CSV export
- Libre CSV export
- Generic CSV (auto-detects columns)

**\*\*Minimum required columns (generic CSV):\*\***

- Timestamp (ISO or epoch minutes)
- Glucose (mg/dL)

**\*\*Optional but recommended:\*\***

- Carbs (grams)
- Insulin (units)
- Notes/events

**\*\*Common import pitfalls:\*\***

- Mixed timezones or missing timezone offsets
- Glucose in mmol/L (convert to mg/dL)



Duplicate timestamps (keep the latest reading)

**\*\*Quick validation:\*\***

Plot glucose vs time to confirm it is smooth and within 40-400 mg/dL

Ensure meal events align with glucose rises (15-90 minutes after)

Nightscout JSON

Dataset registry packs (AIDE, PEDAP, AZT1D, HUPA-UCM)

## 2.9 Generating Custom Reports

```
from iints.analysis.reporting import ClinicalReportGenerator

generator = ClinicalReportGenerator()
generator.generate_pdf(
    results_df, # Your simulation results DataFrame
    safety_report, # Safety summary dict
    "results/custom_report.pdf",
    title="My Algorithm Performance Report",
    include_trend_analysis=True,
    highlight_safety_events=True
)
```

**\*\*Report Contents:\*\***

Glucose trace chart with target range

Insulin delivery chart

Clinical metrics table (TIR, GMI, CV, LBGI, HBGI)

Safety interventions summary

Top intervention reasons

Configuration overview

## 2.10 What Next?

**\*\*Recommended Next Steps:\*\***

[OK] Complete the Getting Started guide

[Metrics] Run benchmark comparisons

[Research] Test with stress scenarios

[Performance] Import real CGM data

[AI] Explore AI predictor integration

[Notes] Generate clinical reports

[Training] Review safety architecture

[Tools] Customize patient profiles

[Package] Package algorithm for distribution

[Announcement] Share results with community

**\*\*Need Help?\*\***

Check Troubleshooting section (9.0)

Review FAQ (17.0)

Join community discussions

[Submit issues on GitHub](#)

## 3. Architecture Overview

### 3.1 System Components

#### IINTS-AF SDK Architecture

Algorithm Layer	Safety Layer	Output Layer
- Custom Algorithm	- Input Validator	- Results
- PID Controller	- Safety Supervisor	- Reports
- ML Predictor	- Safety Config	- Audit
- Pump Emulators		- Metrics

#### Patient Simulation

- Virtual Patient Model
- CGM Sensor Model (noise, lag, dropout)
- Insulin Pump Model (limits, quantization)
- Pharmacokinetics (insulin absorption)
- Pharmacodynamics (glucose response)

### 3.2 Data Flow

1. Algorithm requests insulin dose
2. Input Validator checks glucose values
3. Safety Supervisor applies 9 safety checks
4. Approved dose sent to pump model
5. Pump model simulates delivery (with possible errors)
6. Patient model calculates glucose impact
7. CGM sensor model adds noise/lag
8. New glucose reading returned to algorithm
9. Audit trail logs all decisions
10. Repeat every time step (default: 5 minutes)

### 3.3 Safety Layer Integration

The Independent Safety Supervisor runs **\*\*deterministically\*\*** and can:

- [OK] Override dangerous algorithm requests
- [OK] Log all interventions with reasons
- [OK] Enforce hard limits (hypoglycemia protection)

[OK] Apply dynamic limits (IOB clamping)

[OK] Validate all inputs/outputs

**\*\*Key Principle:\*\*** Safety layer is **\*\*always active\*\*** and cannot be disabled.

## 4. Safety Architecture (Critical Section)

### 4.1 Design Philosophy

**\*\*Safety-First Principles:\*\***

**\*\*Deterministic Overrides\*\***: Same input → same safety decision

**\*\*Fail-Safe Defaults\*\***: When in doubt, reduce insulin

**\*\*Audit Everything\*\***: Every decision logged for accountability

**\*\*Transparent Logic\*\***: Clear reasons for all interventions

**\*\*Configurable Thresholds\*\***: Adapt to different patient profiles

**\*\*Safety Guarantees:\*\***

No algorithm can deliver unsafe doses

Hypoglycemia protection is absolute (< 40 mg/dL emergency stop)

All interventions are logged and explainable

Configuration is validated before simulation starts

### 4.2 SafetyConfig Configuration

```
from iints.core.safety import SafetyConfig

# Default configuration (clinic-safe)
config = SafetyConfig(
    # Hypoglycemia protection
    hypo_cutoff=70.0, # mg/dL - start reducing insulin
    severe_hypo_cutoff=54.0, # mg/dL - emergency stop
    critical_hypo_cutoff=40.0, # mg/dL - immediate termination

    # Hyperglycemia limits
    hyper_cutoff=300.0, # mg/dL - maximum allowed

    # Insulin limits
    max_basal_rate=2.0, # U/hr
    max_bolus=5.0, # U per bolus
    max_iob=10.0, # U total active insulin

    # Rate limits
    max_insulin_per_hour=15.0, # U/hr rolling window
    max_insulin_per_day=80.0, # U/day absolute limit

    # Trend protection
    contract_enabled=True,
    contract_glucose_threshold=90.0, # mg/dL
    contract_trend_threshold=-1.0, # mg/dL per 5 minutes
)
```

**\*\*When to tune SafetyConfig\*\***

Only after you can reproduce a baseline run with stable glucose and no crashes.

Increase strictness (lower cutoffs / lower max rates) when testing new or unstable algorithms.

Relax cutoffs only for controlled research experiments with full audit logs.

**\*\*Recommended baseline ranges (adult research)\*\***

`max\_bolus`: 2-6 U

`max\_basal\_rate`: 1-3 U/hr

`max\_iob`: 6-12 U

`hyper\_cutoff`: 250-300 mg/dL

**\*\*Audit note:\*\*** All SafetyConfig values are written to `run\_metadata.json` and `audit/safety\_summary.json`.

### 4.3 9 Safety Checks Explained

The IndependentSupervisor applies these checks **\*\*in order\*\***:

1. Predictive Hypo Guard [EMERGENCY]
  - If glucose < 70 AND falling fast (-3+ mg/dL per 5min)
  - Action: Suspend insulin for 30 minutes
  - Rationale: Prevent imminent severe hypo
2. Basal Rate Limit [WARNING]
  - If basal > max\_basal\_rate
  - Action: Cap at max\_basal\_rate
  - Rationale: Prevent basal overdose
3. Hard Hypo Cutoff [EMERGENCY]
  - If glucose < 54 mg/dL
  - Action: Suspend all insulin
  - Rationale: Severe hypoglycemia protection
4. Severe Hypo Emergency Stop [EMERGENCY]
  - If glucose < 40 mg/dL
  - Action: Terminate simulation
  - Rationale: Critical hypoglycemia - stop everything
5. Glucose Level Clamp [CRITICAL/WARNING]
  - If glucose > 300 mg/dL
  - Action: Reduce insulin by 50%
  - Rationale: Prevent over-correction
6. Rate-of-Change Trend Stop [CRITICAL]
  - If glucose falling > 5 mg/dL per 5min
  - Action: Suspend insulin
  - Rationale: Rapid drop protection
7. Dynamic IOB Clamp [WARNING]
  - If IOB > max\_iob
  - Action: Reduce dose to stay under max\_iob
  - Rationale: Prevent insulin stacking
8. Bolus Stacking Check [WARNING]
  - If recent boluses > safety limit
  - Action: Delay or reduce bolus
  - Rationale: Prevent bolus overlap

#### 9. 60-Minute Rolling Cap [WARNING]

- If insulin in last 60min > max\_insulin\_per\_hour
- Action: Reduce dose to stay under limit
- Rationale: Hourly limit enforcement

### 4.4 Safety Levels

| Level | Severity | Action |

|-----|-----|-----|

| INFO | Informational | Log only, no intervention |

| WARNING | Potential issue | Adjust dose within safe limits |

| CRITICAL | Serious risk | Significant dose reduction |

| EMERGENCY | Immediate danger | Suspend insulin completely |

### 4.5 Input Validation

The InputValidator checks:

```
# Glucose range validation
if glucose < 20 or glucose > 600:
    raise InvalidGlucoseError(f"Glucose {glucose} outside valid range")

# Insulin request validation
if insulin < 0 or insulin > config.max_bolus:
    raise InvalidInsulinError(f"Insulin {insulin} invalid")

# Timestep validation
if timestep < 1 or timestep > 15:
    raise InvalidTimestepError(f"Timestep {timestep} invalid")
```

### 4.6 Simulation Termination

Automatic termination occurs when:

**\*\*Critical hypoglycemia\*\***: Glucose < 40 mg/dL for 30+ minutes

**\*\*Configuration error\*\***: Invalid safety configuration

**\*\*Algorithm error\*\***: Unhandled exception in algorithm

**\*\*Manual stop\*\***: User interrupts simulation

**\*\*Termination Output\*\***

    `SimulationLimitError` exception raised

    Safety report marks `terminated\_early: true`

    Final glucose and intervention reason logged

    Partial results still saved

## 5. Tutorials and Cookbook

### 5.1 24-Hour Simulation Walkthrough

Complete example from setup to analysis:

```
import iints
import pandas as pd
import matplotlib.pyplot as plt
from iints.core.algorithms.pid_controller import PIDController
```

```

from iints.core.simulator import Simulator, StressEvent

# 1. Setup simulation
sim = Simulator(
    algorithm=PIDController(),
    patient_config="default",
    time_step=5, # 5-minute steps
    enable_profiling=True
)

# 2. Add realistic stress events
sim.add_stress_event(StressEvent(start_time=480, event_type="meal", value=60)) # 8:00 AM
breakfast
sim.add_stress_event(StressEvent(start_time=720, event_type="meal", value=45)) # 12:00 PM
lunch
sim.add_stress_event(StressEvent(start_time=1080, event_type="exercise", value=45)) # 6:00
PM workout
sim.add_stress_event(StressEvent(start_time=1320, event_type="meal", value=75)) # 8:00 PM
dinner

# 3. Run 24-hour simulation
results_df, safety_report = sim.run_batch(duration_minutes=1440)

# 4. Analyze results
print(f"Time in Range (70-180 mg/dL): {iints.metrics.calculate_tir(results_df):.1f}%")
print(f"Glucose Management Indicator: {iints.metrics.calculate_gmi(results_df):.1f}")
print(f"Safety interventions: {safety_report['intervention_count']}")

# 5. Visualize
plt.figure(figsize=(12, 6))
plt.plot(results_df['timestamp'], results_df['glucose_actual_mgdl'])
plt.axhline(180, color='red', linestyle='--', label='Hyperglycemia')
plt.axhline(70, color='green', linestyle='--', label='Target')
plt.axhline(54, color='orange', linestyle='--', label='Hypoglycemia')
plt.title('24-Hour Glucose Profile')
plt.xlabel('Time')
plt.ylabel('Glucose (mg/dL)')
plt.legend()
plt.grid(True)
plt.show()

# 6. Generate report
iints.generate_clinical_report(
    results_df,
    safety_report,
    "results/24hour_report.pdf"
)

```

## 5.2 Building an ML-Hybrid Algorithm

Combine ML prediction with rule-based safety:

```

from iints.api.base_algorithm import InsulinAlgorithm
from iints.research.predictor import load_predictor_service

class MLHybridAlgorithm(InsulinAlgorithm):

```

```

def __init__(self, predictor_path="models/predictor.pt"):
    super().__init__()
    self.predictor = load_predictor_service(predictor_path)
    self.last_prediction = None

def predict_insulin(self, data: AlgorithmInput) -> dict:
    # 1. Get ML prediction (30-min forecast)
    prediction = self.predictor.predict(data)
    self.last_prediction = prediction['glucose_forecast']

    # 2. Rule-based decision with ML insight
    basal = 0.9 # U/hr
    bolus = 0.0

    # 3. Adjust based on prediction
    if prediction['glucose_forecast'] > 200:
        # Aggressive correction if rising
        bolus = (prediction['glucose_forecast'] - 180) / 40
    elif prediction['glucose_forecast'] < 90:
        # Conservative if dropping
        basal = max(0.3, basal * 0.7) # Reduce basal but don't suspend

    return {
        "basal_insulin": basal,
        "bolus_insulin": bolus,
        "reason": f"ML forecast: {prediction['glucose_forecast']:.0f}mg/dL"
    }

```

### 5.3 Running Batch Experiments

Test multiple configurations efficiently:

```

from iints.analysis.batch import run_batch_experiment

configurations = [
    {"algorithm": "PIDController", "patient": "default", "scenario": "baseline"},
    {"algorithm": "PIDController", "patient": "adolescent", "scenario": "meal_challenge"},
    {"algorithm": "MyAlgorithm", "patient": "default", "scenario": "baseline"},
    {"algorithm": "MyAlgorithm", "patient": "adolescent", "scenario": "meal_challenge"},
]

results = run_batch_experiment(
    configurations=configurations,
    duration_minutes=1440,
    output_dir="results/batch_experiment",
    parallel_workers=4 # Use 4 CPU cores
)

# Compare metrics across all runs
comparison_df = results.compare_metrics()
print(comparison_df[['algorithm', 'patient', 'TIR', 'GMI', 'interventions']])

```

### 5.4 Audit Trail Analysis

Every run produces a structured audit trail that explains why the safety layer intervened.

```
import pandas as pd
```

```
audit = pd.read_csv("results/your_run/audit/audit_trail.csv")
print(audit[['timestamp', 'glucose_actual_mgdl', 'action', 'reason']].head())

# Count interventions by type
print(audit['action'].value_counts())
```

#### **\*\*Interpretation tips:\*\***

If `action` is `suspend` or `cap`, the supervisor overrode the algorithm.

If `reason` repeats often, your algorithm is too aggressive for that patient profile.

### 5.5 Custom Safety Thresholds

Use `SafetyConfig` to tighten or relax constraints for research experiments.

```
from iints.core.safety import SafetyConfig
from iints.core.simulator import Simulator

safe_config = SafetyConfig(
    max_bolus=3.0,
    max_iob=8.0,
    hyper_cutoff=250.0,
)

sim = Simulator(algorithm=MyAlgorithm(), safety_config=safe_config)
```

### 5.6 Pump Emulator Benchmarking

Benchmark alternative pump behaviors or commercial-emulator presets.

```
from iints.analysis.hardware_benchmark import benchmark_pump_emulators

bench = benchmark_pump_emulators(duration_minutes=120)
print(bench[['model', 'avg_step_ms', 'max_step_ms']])
```

### 5.7 Live Streaming Simulation

Stream real-time values for dashboards or demos.

```
from iints.core.simulator import Simulator

sim = Simulator(algorithm=MyAlgorithm())
for state in sim.run_stream(duration_minutes=120):
    print(state['timestamp'], state['glucose_actual_mgdl'])
```

### 5.8 Reproducible Runs for Publications

To make results reproducible:

Fix `seed`

Persist `config.json`

Record dataset hash and commit SHA

```
results = iints.run_simulation(
    algorithm=MyAlgorithm(),
    duration_minutes=720,
    seed=123,
)
print(results['run_manifest'])
```



## 6. API Reference

### 6.1 Core Classes

#### InsulinAlgorithm (Abstract Base Class)

```
from iints.api.base_algorithm import InsulinAlgorithm, AlgorithmInput, AlgorithmResult

class MyAlgorithm(InsulinAlgorithm):
    def get_algorithm_metadata(self) -> dict:
        """Return algorithm identification and version info"""
        return {
            "name": "MyAlgorithm",
            "version": "1.0.0",
            "description": "My custom insulin dosing algorithm",
            "author": "Your Name",
            "reference": "Optional citation or paper reference"
        }

    def predict_insulin(self, data: AlgorithmInput) -> dict:
        """
        Calculate insulin dose based on current data

        Args:
            data: AlgorithmInput containing current state

        Returns:
            dict with keys: basal_insulin, bolus_insulin, reason
        """
        # Your algorithm logic here
        return {
            "basal_insulin": 0.9, # U/hr
            "bolus_insulin": 0.0, # U
            "reason": "Stable glucose, maintaining basal rate"
        }

    def reset(self):
        """Reset algorithm state for new simulation"""
        pass
```

#### AlgorithmInput (Dataclass)

```
@dataclass
class AlgorithmInput:
    # Current state
    current_glucose: float # mg/dL
    current_time: datetime # Simulation timestamp

    # Historical context
    glucose_history: List[float] # Last 24 hours (5-min intervals)
    insulin_history: List[float] # Last 24 hours
    carb_history: List[float] # Last 24 hours

    # Calculated values
    iob: float # Insulin on board (U)
    cob: float # Carbs on board (g)
```

```

# Trends
glucose_trend: float # mg/dL per 5 minutes
glucose_acceleration: float # mg/dL per 5 minutes

# Patient info
patient_config: dict # ISF, ICR, basal rates, etc.

# Safety context
last_safety_intervention: Optional[dict] # Last intervention reason

```

## 7. Practical Examples

### 7.1 Complete Algorithm Example

Full working algorithm with comprehensive logic:

```

from iints.api.base_algorithm import InsulinAlgorithm, AlgorithmInput
import numpy as np

class ComprehensiveAlgorithm(InsulinAlgorithm):
    def __init__(self):
        super().__init__()
        self.target_glucose = 120 # mg/dL
        self.isf = 50 # Insulin sensitivity factor
        self.icr = 10 # Insulin-to-carb ratio
        self.basal_rate = 0.9 # U/hr
        self.history = []

    def get_algorithm_metadata(self):
        return {
            "name": "Comprehensive Algorithm",
            "version": "1.0.0",
            "description": "Full-featured algorithm with meal detection and trend analysis"
        }

    def predict_insulin(self, data: AlgorithmInput) -> dict:
        # Store history for trend analysis
        self.history.append(data.current_glucose)
        if len(self.history) > 24:
            self.history.pop(0)

        # Calculate correction bolus
        correction = 0.0
        if data.current_glucose > self.target_glucose:
            correction = (data.current_glucose - self.target_glucose) / self.isf

        # Meal detection (simple version)
        meal_bolus = 0.0
        if data.carb_history and data.carb_history[-1] > 0:
            meal_bolus = data.carb_history[-1] / self.icr

        # Trend adjustment
        trend_adjustment = 0.0
        if data.glucose_trend > 2: # Rising fast
            correction *= 1.2 # More aggressive

```

```

        elif data.glucose_trend < -2: # Dropping fast
            correction *= 0.8 # More conservative

        # IOB safety
        if data.iob > 5: # High IOB
            self.basal_rate = max(0.3, self.basal_rate * 0.8)

        # Final dose calculation
        basal = self.basal_rate
        bolus = correction + meal_bolus

        return {
            "basal_insulin": basal,
            "bolus_insulin": bolus,
            "reason": (
                f"Target {self.target_glucose}mg/dL, "
                f"correction {correction:.2f}U, "
                f"meal {meal_bolus:.2f}U, "
                f"trend {data.glucose_trend:.1f}mg/dL per 5min"
            )
        }

    def reset(self):
        self.history = []

```

## 7.2 CLI vs Python API Comparison

**\*\*Same Task: Run Simulation with Custom Algorithm\*\***

**\*\*CLI Version:\*\***

```

# Create algorithm
iints new-algo --name MyAlgorithm --output-dir algorithms/

# Edit algorithm file
nano algorithms/my_algorithm.py

# Run simulation
iints run \
    --algo algorithms/my_algorithm.py \
    --patient-config-name default \
    --scenario-path scenarios/meal_challenge.json \
    --duration 1440 \
    --output-dir results/cli_run

```

**\*\*Python Version:\*\***

```

from iints import run_simulation
from iints.core.algorithms.custom import MyAlgorithm

results = run_simulation(
    algorithm=MyAlgorithm(),
    patient_config="default",
    scenario="scenarios/meal_challenge.json",
    duration_minutes=1440,
    output_dir="results/python_run"
)

```

**\*\*When to Use CLI:\*\***

- Quick testing and iteration
- Batch processing multiple scenarios
- Integration with shell scripts
- CI/CD pipelines

**\*\*When to Use Python API:\*\***

- Complex algorithm development
- Integration with other Python tools
- Custom analysis pipelines
- Jupyter notebook exploration

### 7.3 Stress Testing Patterns

**\*\*Pattern 1: Meal Challenge Test\*\***

```
# Test algorithm response to large meal
sim.add_stress_event(StressEvent(
    start_time=480, # 8:00 AM
    event_type="meal",
    value=100 # 100g carbs (large pizza meal)
))
```

**\*\*Pattern 2: Exercise Stress Test\*\***

```
# Test algorithm response to exercise
sim.add_stress_event(StressEvent(
    start_time=1080, # 6:00 PM
    event_type="exercise",
    value=60 # 60 minutes intense exercise
))
```

**\*\*Pattern 3: Sensor Noise Test\*\***

```
# Test algorithm robustness to CGM noise
sim.add_stress_event(StressEvent(
    start_time=300, # 5:00 AM
    event_type="sensor_noise",
    value=20.0 # 20 mg/dL standard deviation
))
```

**\*\*Pattern 4: Combined Stress Test\*\***

```
# Realistic day with multiple stressors
sim.add_stress_event(StressEvent(420, "meal", 45)) # 7:00 AM breakfast
sim.add_stress_event(StressEvent(660, "meal", 60)) # 11:00 AM snack
sim.add_stress_event(StressEvent(900, "exercise", 30)) # 3:00 PM walk
sim.add_stress_event(StressEvent(1020, "meal", 80)) # 5:00 PM dinner
sim.add_stress_event(StressEvent(1200, "sensor_noise", 15.0)) # 8:00 PM sensor issues
```

### 7.4 Human-in-the-Loop Integration

Add manual interventions during simulation:

```
def human_in_loop_callback(context):
    """
    Called at each simulation step.
    Return dict with manual interventions or None.
    """
```

```

# Example: Rescue carbs for hypoglycemia
if context["glucose_actual_mgdl"] < 65:
    return {
        "additional_carbs": 15, # 15g fast-acting carbs
        "note": "Human intervention: rescue carbs for hypo"
    }

# Example: Manual bolus correction
if context["glucose_actual_mgdl"] > 250:
    return {
        "additional_bolus": 1.5, # 1.5U correction
        "note": "Human intervention: manual correction bolus"
    }

# Example: Suspend insulin before exercise
if context["time"] > 1080 and context["time"] < 1140: # 6-7 PM
    return {
        "suspend_insulin": True,
        "note": "Human intervention: exercise suspension"
    }

return None # No intervention

# Use in simulator
sim = Simulator(
    algorithm=MyAlgorithm(),
    patient_config="default",
    on_step=human_in_loop_callback
)

```

## 7.5 Data Import Workflows

### **\*\*Workflow 1: Dexcom CSV Import\*\***

```

# CLI import
iints import-data \
  --input-csv dexcom_export.csv \
  --data-format dexcom \
  --output-dir results/dexcom_import

# Use imported scenario
iints run \
  --algo algorithms/my_algorithm.py \
  --scenario-path results/dexcom_import/scenario.json

```

### **\*\*Workflow 2: Nightscout Import\*\***

```

# Install Nightscout extra
pip install iints-sdk-python35[nightscout]

# Import from Nightscout
iints import-nightscout \
  --url https://your-nightscout-site.herokuapp.com \
  --days 7 \
  --output-dir results/nightscout_import

```

### **\*\*Workflow 3: Dataset Registry\*\***

```

# List available datasets

```

```
iints data list

# Fetch specific dataset
iints data fetch aide_t1d --output-dir data_packs/aide

# Use in simulation
iints run \
  --algo algorithms/my_algorithm.py \
  --scenario-path data_packs/aide/scenarios/patient_001.json
```

## 8. Advanced Topics

### 8.1 Commercial Pump Emulators

Test against real pump behavior:

```
from iints.emulation import Medtronic780G, Omnipod5, TandemControlIQ

# Compare your algorithm against commercial pumps
pumps = [
    ("MyAlgorithm", MyAlgorithm()),
    ("Medtronic 780G", Medtronic780G()),
    ("Omnipod 5", Omnipod5()),
    ("Tandem Control-IQ", TandemControlIQ())
]

results = {}
for name, pump_algo in pumps:
    results[name] = run_simulation(
        algorithm=pump_algo,
        patient_config="default",
        duration_minutes=1440
    )

# Compare metrics
compare_metrics(results)
```

### 8.2 Dataset Registry Usage

Access real-world datasets:

```
# List all available datasets
iints data list

# Get info about specific dataset
iints data info aide_t1d

# Fetch dataset
iints data fetch aide_t1d --output-dir data_packs/aide

# Cite dataset in publication
iints data cite aide_t1d
```

**\*\*Available Datasets:\*\***

```
`aide_t1d`: AIDE Type 1 Diabetes Dataset
`azt1d`: Arizona Type 1 Diabetes Dataset
`ohio_t1dm`: OhioT1DM Dataset
```

`sample`: Bundled demo data (no download needed)

### **\*\*Integrity and reproducibility\*\***

Every dataset entry includes a SHA-256 checksum and citation metadata.

The fetch command validates the checksum automatically.

Use `iints data info <dataset>` to record version and hash in your paper.

### **\*\*Typical layout after fetch\*\***

```
data_packs/
  public/
    ohio_t1dm/
      raw/...
      processed/...
```

## 8.3 Reproducibility Techniques

Ensure identical results across runs:

```
# Method 1: Set random seed
results = run_simulation(
    algorithm=MyAlgorithm(),
    seed=42, # Fixed random seed
    patient_config="default"
)

# Method 2: Use deterministic patient profile
profile = PatientProfile(
    isf=50,
    icr=10,
    basal_rate=0.9,
    dawn_phenomenon=0.3, # Fixed dawn effect
    seed=42
)

# Method 3: SHA-256 verification
from iints.utils.hashing import verify_manifest

if verify_manifest("results/run_001/run_manifest.json"):
    print("Results verified - not tampered with")
```

## 8.4 Performance Profiling

Measure algorithm performance:

```
sim = Simulator(
    algorithm=MyAlgorithm(),
    patient_config="default",
    enable_profiling=True
)

results_df, safety_report = sim.run_batch(duration_minutes=1440)

# Access performance data
profiling = safety_report["performance_report"]
print(f"Algorithm latency: {profiling['algorithm_latency_ms']:.2f}ms")
print(f"Safety latency: {profiling['safety_latency_ms']:.2f}ms")
print(f"Total steps: {profiling['total_steps']}")
```

```
print(f"Total time: {profiling['total_time_s']:.2f}s")
```

## 8.5 Custom Metrics Calculation

Define your own performance metrics:

```
def calculate_custom_metric(results_df):
    """Calculate custom performance metric"""

    # Time in tight range (80-140 mg/dL)
    tight_range = ((results_df['glucose_actual_mgdl'] >= 80) &
                   (results_df['glucose_actual_mgdl'] <= 140)).mean() * 100

    # Glucose variability score
    cv = results_df['glucose_actual_mgdl'].std() /
        results_df['glucose_actual_mgdl'].mean() * 100

    # Hypoglycemia risk score
    hypo_risk = (results_df['glucose_actual_mgdl'] < 70).sum() / len(results_df)

    # Hyperglycemia risk score
    hyper_risk = (results_df['glucose_actual_mgdl'] > 250).sum() / len(results_df)

    # Composite score (lower is better)
    composite_score = (100 - tight_range) + cv + (hypo_risk * 100) + (hyper_risk * 50)

    return {
        'tight_range_percent': tight_range,
        'cv_percent': cv,
        'hypo_risk_percent': hypo_risk * 100,
        'hyper_risk_percent': hyper_risk * 100,
        'composite_score': composite_score
    }

# Use with your results
metrics = calculate_custom_metric(results_df)
print(f"Custom Score: {metrics['composite_score']:.1f}")
```

## 9. Troubleshooting

### 9.1 Installation Issues

**\*\*Issue: ModuleNotFoundError after installation\*\***

**\*\*Solution:\*\***

```
# Make sure you're using the correct Python environment
which python # Should show your virtual environment path
```

```
# Reinstall in development mode if needed
pip install -e .
```

**\*\*Issue: pip install fails with dependency errors\*\***

**\*\*Solution:\*\***

```
# Upgrade pip first
pip install --upgrade pip setuptools wheel
```

```
# Try installing with --no-cache-dir
```



```
pip install --no-cache-dir iints-sdk-python35
```

```
# For specific errors, check Python version
python3 --version # Must be 3.10+
```

## 9.2 Simulation Issues

**\*\*Issue: Simulation runs very slowly\*\***

**\*\*Solution:\*\***

```
# Increase time step (default is 5 minutes)
sim = Simulator(time_step=10) # 10-minute steps

# Disable profiling if not needed
sim = Simulator(enable_profiling=False)

# Reduce duration for testing
results = sim.run_batch(duration_minutes=720) # 12 hours instead of 24
```

**\*\*Issue: Simulation terminates early\*\***

**\*\*Solution:\*\***

```
# Check safety report for termination reason
print(safety_report['termination_reason'])

# Common causes:
# - Critical hypoglycemia (< 40 mg/dL for 30+ minutes)
# - Algorithm exception
# - Configuration error

# Adjust safety limits if needed
from iints.core.safety import SafetyConfig
config = SafetyConfig(critical_hypo_cutoff=35.0) # Lower threshold
```

## 9.3 Algorithm Development Issues

**\*\*Issue: Algorithm not appearing in CLI\*\***

**\*\*Solution:\*\***

```
# Make sure algorithm is properly registered
# 1. Inherits from InsulinAlgorithm
# 2. Implements all required methods
# 3. Has valid get_algorithm_metadata()

# Check with:
from iints.api.registry import list_algorithms
print(list_algorithms())
```

**\*\*Issue: Safety supervisor blocking all doses\*\***

**\*\*Solution:\*\***

```
# Check safety configuration
print(safety_report['config'])

# Common issues:
# - max_basal_rate too low
# - hypo_cutoff too high
# - contract enabled with aggressive settings
```

```
# Adjust configuration:
config = SafetyConfig(
    hypo_cutoff=60.0, # Higher threshold
    max_basal_rate=1.5 # Higher limit
)
```

## 9.4 Data Import Issues

**\*\*Issue: CSV import fails\*\***

**\*\*Solution:\*\***

```
# Check CSV format
# Required columns: timestamp, glucose
# Optional: carbs, insulin

# Try auto-detect format
data = import_cgm_csv("your_file.csv", data_format="auto")

# Specify column names manually if needed
data = import_cgm_csv(
    "your_file.csv",
    data_format="generic",
    timestamp_col="Date",
    glucose_col="BG",
    carbs_col="Carbs"
)
```

## 9.5 Performance Issues

**\*\*Issue: High memory usage\*\***

**\*\*Solution:\*\***

```
# Reduce history size
sim = Simulator(max_history_hours=12) # Default is 24

# Disable audit trail if not needed
sim = Simulator(enable_audit=False)

# Process in batches
for i in range(10):
    results = sim.run_batch(duration_minutes=144) # 2.4 hours per batch
    process_results(results)
```

# 10. Quick Reference

## 10.1 Essential CLI Commands

```
# Initialize project
iints init --project-name my_project

# Quickstart
iints quickstart --project-name demo

# Run simulation
iints run --algo algorithms/my_algo.py --duration 1440
```

```
# Run with preset
iints presets run --name baseline_t1d --algo algorithms/my_algo.py

# Benchmark
iints benchmark --algo-to-benchmark algorithms/my_algo.py

# Import data
iints import-data --input-csv my_cgm.csv --data-format dexcom

# List presets
iints presets list

# Generate scenario
iints scenarios generate --name "Stress Test"

# Validate files
iints validate --scenario-path scenarios/my_scenario.json
```

## 10.2 Python Code Snippets

### **\*\*Minimal Simulation:\*\***

```
from iints import run_simulation
from iints.core.algorithms.pid_controller import PIDController

results = run_simulation(
    algorithm=PIDController(),
    duration_minutes=720
)
```

### **\*\*Custom Algorithm:\*\***

```
from iints.api.base_algorithm import InsulinAlgorithm

class MyAlgorithm(InsulinAlgorithm):
    def predict_insulin(self, data):
        return {"basal_insulin": 0.9, "bolus_insulin": 0.0}
```

### **\*\*Load Results:\*\***

```
import pandas as pd
df = pd.read_csv(results['results_csv'])
print(df[['timestamp', 'glucose_actual_mgdl', 'insulin_delivered']].head())
```

## 10.3 Safety Thresholds (Defaults)

```
SafetyConfig(
    hypo_cutoff=70.0,          # Start reducing insulin
    severe_hypo_cutoff=54.0,   # Emergency suspension
    critical_hypo_cutoff=40.0, # Terminate simulation
    hyper_cutoff=300.0,        # Maximum glucose
    max_basal_rate=2.0,        # U/hr
    max_bolus=5.0,             # U per bolus
    max_iob=10.0,              # U total active insulin
    max_insulin_per_hour=15.0, # U/hr rolling window
    max_insulin_per_day=80.0   # U/day absolute limit
)
```

## 10.4 Clinical Metrics Targets (ATTD/ADA)

Metric	Target Range
TIR (70-180 mg/dL)	>70%
TBR (<70 mg/dL)	<4%
TBR (<54 mg/dL)	<1%
GMI	<7.0%
CV	<36%
LBGI	<1.1
HBGI	<5.0

## 11. Glossary

**\*\*Algorithm\*\***: Code that calculates insulin doses based on current and historical data

**\*\*Basal Insulin\*\***: Background insulin delivery rate (U/hr)

**\*\*Bolus Insulin\*\***: Additional insulin for meals or corrections (U)

**\*\*CGM\*\***: Continuous Glucose Monitor - device that measures glucose every 5 minutes

**\*\*COB\*\***: Carbs On Board - carbohydrates still being absorbed

**\*\*GMI\*\***: Glucose Management Indicator - estimate of A1C from CGM data

**\*\*IOB\*\***: Insulin On Board - insulin still active in the body

**\*\*ISF\*\***: Insulin Sensitivity Factor - how much 1U of insulin lowers glucose (mg/dL)

**\*\*ICR\*\***: Insulin-to-Carb Ratio - grams of carbs covered by 1U of insulin

**\*\*TIR\*\***: Time In Range - percentage of time in target glucose range (70-180 mg/dL)

**\*\*TBRL1\*\***: Time Below Range Level 1 - percentage of time <70 mg/dL

**\*\*TBRL2\*\***: Time Below Range Level 2 - percentage of time <54 mg/dL

**\*\*TAR\*\***: Time Above Range - percentage of time >180 mg/dL

**\*\*CV\*\***: Coefficient of Variation - measure of glucose variability

**\*\*LBGI\*\***: Low Blood Glucose Index - measure of hypoglycemia risk

**\*\*HBGI\*\***: High Blood Glucose Index - measure of hyperglycemia risk

## Need More Help?

### **\*\*Documentation:\*\***

Comprehensive Guide: ``docs/COMPREHENSIVE_GUIDE.md``

Technical README: ``docs/TECHNICAL_README.md``

API Stability: ``API_STABILITY.md``

### **\*\*Community:\*\***

GitHub Issues: <https://github.com/python35/IINTS-SDK/issues>

Discussion Forum: [Link]

Email Support: [support@iints.org](mailto:support@iints.org)

### **\*\*Citing IINTS-AF:\*\***

```
@software{IINTS_AF_SDK,
  author = {Bobbbaers, Rune},
  title = {IINTS-AF: Intelligent Insulin Titration System for Artificial Pancreas},
  year = {2026},
```

```
version = {0.1.19},  
url = {https://github.com/python35/IINTS-SDK},  
note = {Pre-clinical research software for insulin dosing algorithm validation}  
}
```

**\*\*PRE-CLINICAL USE ONLY - NOT FOR PATIENT CARE\*\***

This document is provided for research and educational purposes only.

The IINTS-AF SDK is not a medical device and should not be used for clinical decision-making.

© 2026 IINTS-AF Project. Licensed under MIT License.